

Product Name :
TVET Press Laminate and Wood Coating

Product Code :
NLAB-EDUCATIONTVET20005



Description :

TVET Press Laminate and Wood Coating

Technical Specification :

Max nipping area: 430x320mm/ larger than A3 size

Max press weight(force): 1500kg

Max wood thickness: 13cm

Machine size: 36x46x50cm

Machine weight: 25kg

Naugalabequipments

Website: www.naugalabequipments.com, **Email:** sales@naugalabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003