

Product Name :
TVET Sanitation Fittings Training Panel

Product Code :
NLAB-EDUCATIONTVET27001



Description :

TVET Sanitation Fittings Training Panel

Technical Specification :

investigating the most common sanitation fittings

facilities for hot and cold water supply using quick-release couplings

possibility to test different fittings simultaneously

flushing valve

two handle mixer for a shower

two handle mixer as water tap

dial thermometers

Bourdon tube manometers

rotameters

adjustable pressure reducing valves

Technical data

Flushing valve

1,2...5bar

Pressure reducing valve

max. initial pressure: 25bar

outlSR pressure: 1,5...6bar, up to 70°C

Measuring ranges

temperature: 2x 0...60°C, 1x 0...100°C

pressure: 2x 0...6bar

flow rate: 160...1600L/h

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003