

**Product Name :**  
TVET Soil Water Potential Locator

**Product Code :**  
NLAB-EDUCATIONTVET41009



**Description :**

TVET Soil Water Potential Locator

**Technical Specification :**

High accuracy and high resolution

Easy to operate, portable, flexible, able to purchase one set with multi-probes.

Automatically get the peak data of water potential

Time, peak data, automatic record function.

Record temperature, water potential, time and stored No at the same time.

Automatically turns off function

Max pressure: 100kpa

Resolution: 0.01kpa

Accuracy:  $\pm 1\%$  FS

Temperature testing range:  $-55+150^{\circ}\text{C}$

Temperature accuracy:  $\pm 0.5^{\circ}\text{C}$

Qualities of standard probe: 3

soil temperature index

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003