

**Product Name :**  
TVET Sorting and Packaging Machine Trainer

**Product Code :**  
NLAB-EDUCATIONTVET29002



**Description :**

TVET Sorting and Packaging Machine Trainer

**Technical Specification :**

**Brief Introduction**

automatic granule packing machine is mainly used packing different kinds of coarse powder?granule material side sealing flat bag or pillow bag filling and packing. Specific materials that can be packed: all kinds of soup powder, seasonng, coffee, soybean, solt, melon seeds, tea, etc..

**Products Composition and Packing Process**

**General Composition**

The machine is composed by feeding barrel, rotary measuring cup feeding equipment, thin film supply equipment, thin film molding machine, vertical heat sealing equipment, horizontal heat sealing cutting equipment, electric-control part, machine body part equipment etc.. All the parts that contact with packing material is stainless steel.

**Packing technological process**

Putting Film & Bag Forming?heat Sealing?Filling?Bag Cutting?

Finishing Packing

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003