

Product Name :
TVET Sousvide

Product Code :
NLAB-EDUCATIONTVET22004



Description :

TVET Sousvide

Technical Specification :

Power: Approx. 1200 Watt

Temperature range: Approx. 0-95°C

Dimensions: Min. 37 x 18 x 7 cm

Weight: Approx. 1 kg

AC: 220-230V, 60Hz

Safety: Thermal fuse

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

