

Product Name :
TVET Speed and Position Transducer and Control Trainer

Product Code :
NLAB-EDUCATIONTVET25001



Description :

TVET Speed and Position Transducer and Control Trainer

Technical Specification :

DC permanent magnet motor; general notation and mathematical moderlization; Electriocal and mechanical chrastarestics

Detection of chrastarestic of a mature reaction DC m otor

Detection cof chrastaretic optoelectronic speed tradecers

Transducer; potentiometric transducers;

The study of response of proprtlonal?intergrative and derivertive controllers at diffrent input signals

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003