

Product Name :
TVET Standard Triaxial Testing Apparatus

Product Code :
NLAB-EDUCATIONTVET40014

Description :

TVET Standard Triaxial Testing Apparatus

Technical Specification :

Speciament size 39.1mm,Ø61.8mm,Ø101mm

Capacity 60kN

LCD digital display

confining pressure 0Mpa-2Mpa

back pressure 0Mpa-2Mpa

Deformation volume 0ml-50ml, 0ml- 100ml

Working table travel: minimum 0- 100 mm

Axial force 0 - 10kN,0 - 30kN,0 - 60kN

Pore ??pressure 0mPa-2MPa

Sear speed 0.0001mm /min-2.4mm/min.0.0001mm/min-4.8mm/min

Power: approx 220V, 50Hz

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003