

**Product Name :**  
TVET Steam Ironing Table

**Product Code :**  
NLAB-EDUCATIONTVET21008



**Description :**

TVET Steam Ironing Table

**Technical Specification :**

Iron table:

With function of vacuum

Voltage: 220V

Motor power: approx. 550 W

Max. Pressure of air: approx. 150 Pa

Include ironer

Steam generators:

Voltage: 220V

Electricity heating power: approx. 3 kW

Rated working pressure: approx. 0.4 MPa

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003