

**Product Name :**  
TVET Total Total Station

**Product Code :**  
NLAB-EDUCATIONTVET40007



**Description :**

TVET Total Total Station

**Technical Specification :**

Angle Measurement(Hz,V)

Measurement Method Absolute Encoding

Minimum Readout 1"/5"/10"

Angle Accuracy 2?

Distance measurement (with reflector)

Single prism 2500m

Reflective Sheet 800m

Accuracy 2mm+2ppm

Measuring time(Fine/Quick/Tracking) 1.5s/1s/0.5s

Distance measurement (reflectorless)

Reflectorless Range 350m

Single Prism 7500m

Accuracy 3mm+2ppm

Measuring time 1.5s

Telescope

Magnification 30X

Field of View 1°30'

Minimum Focus Distance 1.5m

Illuminated reticle

Compensator

System Single-axis liquid tilt sensor / Dual-axis (optional)

Working Range  $\pm 3'$

Setting Accuracy 1'

Communication

Interface Standard RS232, SD card, USB drive, mini-usb

Internal Data Memory Approx. 20,000 points

Operation system Real-time operation system

Display : High resolution backlight black and white display with contrast adjustable/Graphics:192x96 pixels,  
Character:6 lines x 25 characters

Keyboard 2 sides Alphanumeric backlit crystal keyboard

Standard Accessories:

Tripod, Instruction Manual, Reflector

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003