

**Product Name :**  
TVET Vapor Pressure Tester

**Product Code :**  
NLAB-EDUCATIONTVET39004



**Description :**

TVET Vapor Pressure Tester

**Technical Specification :**

Floor type structure. Stainless steel working table. It is good looking and easy to operate. Accessories are complete.

Digital temperature controller. The temperature in the water bath is uniform. Temperature control accuracy is high. It can meet the requirements of test.

Two precise pressure meters to detect the gas pressure inside the bomb.

Two sets of vapor pressure bombs. It can do synchronous tests.

Power supply: AC(220±10%)V, 50Hz

Heating power of bath: 1600W

Temperature control range of bath: Ambient to 90°C

Temperature control accuracy of bath: ±0.1°C

Accuracy of pressure meter: ±0.4%

Ambient temperature: (-10~35)°C

Relative humidity: 85%

Maximum power consumption: 1700W

Dimension: 350mm×340mm×750mm

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003