

Product Name :
Twin Rack

Product Code :
NLAB-EDUSCHOOLEQ135002



Description :

Twin Rack

Technical Specification :

Twin rack for both 0.5 ml. & 1.5 ml. micro centrifuge tubes.

Each side has 108 holes for 0.5 & 1.5 ml. MCT's respectively.

The transparent top provides an excellent view as well as a protection cover for the tubes stored.

Each face of the base has a unique alpha-numeric index for easy identification of the tube.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003