

Product Name :
Two Door Vertical Freezer

Product Code :
NLAB-VOCATIONALTJET27053

Description :

Two Door Vertical Freezer

Technical Specification :

Body Material: Stainless Steel

Capacity: 500 Litres

Defrost Type: Manual Defrost

Power source: Electric

Size: Large

Voltage: 310 V

Frequency: 60 Hz

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003