

Product Name :
Two Gang Prism Mould 25x25x285 mm to ASTM

Product Code :
NLAB-CEMETLAB130021



Description :

Two Gang Prism Mould 25x25x285 mm to ASTM

Technical Specification :

Two gang prism mould 25x25x285 mm complete with contact points. Gauge length 250mm, total length 285 mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003