

Product Name : Two Tank Fryer	Product Code : NLAB-VOCATIONALTVET27054
Description :	
Two Tank Fryer	
Technical Specification :	
Material: Stainless Steel	
Dimension: 580x440x310 mm	
Color: Silver	
Manual Thermostat Temperature: 0-200 Degree Celsius	
Features:	

Reliable

Optimum quality

Perfect finish

Precise design

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003