

Product Name :
U Shaped Alnico Magnet

Product Code :
NLAB-ENGINEERINGLB12003



Description :

U Shaped Alnico Magnet

Technical Specification :

U Shaped Alnico Magnet Lab Equipments Manufacturer, Suppliers & Exporters

Half red & half blue (for N & S) spray painted, with keepers. Made of cobalt steel 12 x 6 mm - available in different sizes.

Bulk Manufacturers & Suppliers of U Shape Magnet Alnico with High Power. Exporters of U Shaped Alnico Magnets

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003