

Product Name :
Ultrasonic Scaler

Product Code :
NLABMEDICAL-M32100016

Description :

Ultrasonic Scaler India, Ultrasonic Scaler manufacturer India, suppliers, Medical Ultrasonic Scaler India, Hospital Ultrasonic Scaler India

Technical Specification :

Ultrasonic Scaler is a powerful piezoelectric scaler with stable ultrasonic oscillation of 32KHz with very wide energy scope

Titanium Hand piece is autoclavable

Great patient comfort with high frequency and low amplitude

Hand piece: never heat up after long time use

Efficient performance to remove intractable calculus

Fine tuning control on power supply and water spray

Endodontic Function:

Turn the power below half of Top Selector and replace tip with endo kit

(Endoholder = U-file) to operate unit for root canal preparation.

Smoother canal walls grant better obturation

Rapid and superior operating for root canal treatment

Root canal enlargement and swirling irrigation for debridement

Minimal dentin removal keeps structural integrity of the tooth

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003