

Product Name :
Ultrasonic Screen Cleaning Device

Product Code :
NLAB-EDUCATIONTVET13016



Description :

Ultrasonic Screen Cleaning Device

Technical Specification :

Analog setting

Stainless steel tank

Ultrasound peak power: 480 W

SweepTec modulation

Timer: 1 to 15 minutes and continuous walking

Frequency: 35 kHz, automatic adjustment

Dimensions of the tank: Ø24 x 13 cm

Capacity: 5.6 liters (necessary working filling 4.0 liters)

Power supply: 220V - 50Hz

Insertion basket, type K 6, stainless steel, Ø21.6 x 5 cm, mesh 5 x 5 mm

Cover, type D 6

Bath liquid, 25 liters

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003