

**Product Name :**  
Ultratek Full Body Harness

**Product Code :**  
NLAB-SAFETYLABP20034



**Description :**

Ultratek Full Body Harness

**Technical Specification :**

Special Features: Multipurpose harness, ideal for professional use and rescue One external attachment point for rescue & descending Ergonomic Horizontal H-shaped thigh straps

Facilitate vertical and horizontal movements.

Breathable open-core padding for extra comfort in the back and thigh Automatic buckle for quick adjustment Easy no-tangle donning

No discomfort or pressure on the groin.

Improved comfort in case of a fall.

Anchor Point: One dorsal D-ring, Two chest attachment textile loop, One ventral loop, Two lateral D-ring.

Materials: Fray-proof dope-dyed polyester webbing Metal quick buckles of high strength alloy steel.

Stitching: High-strength polyester computer-controlled stitching ensures consistency.

Adjustability: Size adjustment straps at shoulder, waist and thighs.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003