

Product Name :
Ultratek Full Body Harness

Product Code :
NLAB-SAFETYLABP20034



Description :

Ultratek Full Body Harness

Technical Specification :

Special Features: Multipurpose harness, ideal for professional use and rescue One external attachment point for rescue & descending Ergonomic Horizontal H-shaped thigh straps

Facilitate vertical and horizontal movements.

Breathable open-core padding for extra comfort in the back and thigh Automatic buckle for quick adjustment
Easy no-tangle donning

No discomfort or pressure on the groin.

Improved comfort in case of a fall.

Anchor Point: One dorsal D-ring, Two chest attachment textile loop, One ventral loop, Two lateral D-ring.

Materials: Fray-proof dope-dyed polyester webbing Metal quick buckles of high strength alloy steel.

Stitching: High-strength polyester computer-controlled stitching ensures consistency.

Adjustability: Size adjustment straps at shoulder, waist and thighs.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003