

**Product Name :**  
Umbilical Cord Clamp

**Product Code :**  
NLABMEDICAL-M3020001



**Description :**

Umbilical Cord Clamp

**Technical Specification :**

Umbilical Cord Clamp is used for clamping the umbilical cord of newly born baby.

Designed for clamping the umbilical cord immediately after the birth.

Manufactured from non toxic, medical grade polymer.

All grooves are design with the length.

Prevents accidental re-opening.

Finger grip ensures safe & convenient handling.

Provides suitable direction for clamping with smooth re-open facilities.

Sterile, individually packed in peelable pouch pack.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003