

**Product Name :**  
Unconfined Compression Test Apparatus

**Product Code :**  
NLAB-COMPSLAB280002



**Description :**

Unconfined Compression Test Apparatus

**Technical Specification :**

Dimensions: 38 mm x 76 mm

Automation Grade: Automatic

Finishing: Galvanized

Frequency (Hertz): 50

No of Load springs: 4

Load Springs Weight: 0.2 kg, 0.4 kg, 0.8 kg and 1.6 kg/mm

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

