

**Product Name :**  
Unconfined Compression Tester

**Product Code :**  
NLAB-EDUCATIONTVET40006



**Description :**

Unconfined Compression Tester

**Technical Specification :**

Max.capacity: 0.6KN

Dimension of specimen: 39.1mm x 80mm

Speed ??: 2.4mm/min

Dimensions 500 x 170 x 380 mm(LWH)

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003