

Product Name :
Unconfined Compression Testing Apparatus

Product Code :
NLAB-SOILELAB140011



Description :

Unconfined Compression Testing Apparatus

Technical Specification :

Comprises a screw operated load frame, cap.5000kg with a gear box and motor drive giving 1.25, 1.5, 2.5mm/min. Rates of strain, a pair of cone seating, adaptor for proving ring, and stain dial gauge bracket. Supplied with one pair of male/female coning tools for 38mm dia samples but without ring and dial gauge. Suitable for operation on 230 V, Single Phase A.C. Supply.

Accessories:

Coning tools in pairs (Male & Female) for samples having diameter 38, 50, 75 & 100mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003