

Product Name : Undercounter Blood Bank Refrigerators	Product Code : NLABMEDICAL-M3200010
--	---

--

Description :
Undercounter Blood Bank Refrigerators India, Undercounter Blood Bank Refrigerators manufacturer India, suppliers, Medical Undercounter Blood Bank Refrigerators India, Hospital Undercounter Blood Bank Refrigerators India

Technical Specification :
Factory set to +3.5°C (+38°F).
Maintains +1° / -1° C product temperature.
Heavy-duty 1/4 hp compressor.
Off-cycle automatic defrost.
Standard with temperature surveillance module (includes alarm, digital temperature display and 6 chart recorder).
Two durable, stainless steel "floating" drawers storing up to (50) 450 ml bags.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003