

**Product Name :**  
Universal Dust Collector

**Product Code :**  
NLAB-MACHINELABWK33028



**Description :**

Universal Dust Collector

**Technical Specification :**

Chip and saw dust extraction

Motor Power: 1900W

Air flow: 1600 m3/hr

Hose Diameter: 1 x 150mm

Bag capacity: 150 litres

Come complete with 785 x 940mm Waste Sack Orange (Pkt 10);

Micron filter – 30 micron

Cartridge Filter

Wood Workshop Portable Tools Universal Dust Collector Manufacturer, Supplier and Universal Dust Collector Exporter in India

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003