

Product Name :
Universal Extrudes (Hydraulic)

Product Code :
NLAB-TESTINGLABM109037



Description :

Universal Extrudes (Hydraulic)

Technical Specification :

The extractor frame is used for taking out soil samples compacted or undisturbed, from 100mm dia and 150mm dia cylinders such as Core Cutters, Proctor moulds, C.B.R. moulds etc.

Specification:

It consists of a 50Kn capacity hand operated, hydraulic jack mounted on a suitable frame. Two plungers, one for 4"/100mm dia and the other of 6"/150mm dia moulds are supplied Height of thrust plate is adjustable. Set of plungers adaptors and thrust plates for 38mm, 50mm and 75mm dia specimen.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003