

**Product Name :**  
Universal Milling Machine

**Product Code :**  
NLAB-MACHINELABWK13002



**Description :**

Universal Milling Machine

**Technical Specification :**

**TABLE**

**1250**  
Surface

**3x13**  
Number & Size  
of T-Slots

**50**  
Distance  
between T-  
Slots

## FEEDS

~~300~~ss

~~300~~tical

~~300~~gitudinal

## SPINDLE

~~250~~r

~~0(50 to 360)~~  
Speeds

~~ISO 40~~

~~10-300~~ce from  
Spindle to Top  
(min. - max.)

## MOTORS

~~3HP~~le  
1440RPM

~~0.10~~HP. 2800  
RPM

**WEIGHT**

~~2500 Kg.~~

~~2500 Kg.~~

**VERTICAL AT  
TACHMENT  
ON RAM**

~~3500~~

~~7000~~

**6** speeds for  
Vertical Head  
~~(9000)~~e

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003