

Product Name :
Universal Moisture Tester

Product Code :
NLAB-SOILELAB140034



Description :

Universal Moisture Tester

Technical Specification :

Usage/Application: Laboratory

Type: Digital Grain Moisture Meter

Measuring Range: 0 to 70%

Temperature Range: 10 to 60 Degree Celsius

Material: Mild Steel

Display Type: Digital

Riding on unfathomable volumes of industrial expertise, we are providing a broad array of Universal Digital Moisture Meter.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments