

Product Name :
Universal Oven, Type UN260

Product Code :
NLAB-EDUCATIONTVET15022



Description :

Universal Oven, Type UN260

Technical Specification :

Ventilation by natural convection

Stainless steel construction

Temperature: (ambient + 5 ° C) to 300 ° C

Capacity: 256 liters

PID regulation, 1 minute timer to 99 days and 23 hours

High resolution TFT color display

Slides for tray / grid: 9

Internal dimensions (WxHxD): 640 x 800 x 500 mm

External dimensions (WxHxD): 824 x 1,186 x 684 mm

Power supply: 230V - 50 / 60Hz - 3.400W

Weight: 110 kg

Supplied with 2 stainless steel grids and factory calibration certificate at 160 ° C

Additional grid for model 260

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003