

<b>Product Name :</b> Universal Stereoscopic Binocular Microscope	<b>Product Code :</b> NLAB-EDUSCHOOLEQ1730139
<b>Description :</b>	
Universal Stereoscopic Binocular Microscope	
<b>Technical Specification :</b>	
Eyepieces Standard extra wide field eyepieces WF-10X is provided.	
Focusing By Rack and Pinion mechanism with adjustable tension collar.	
Magnification Standard magnification 20X to 60X Base Multi-purpose Universal Base Stand Objectives 2 sets of objectives 2X and 4X, selected by rotating nose piece through 90-degrees.	
Observation Head Binocular head inclined at 45-deg., rotated through 360-deg with adjustable inter-papillary distance.	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003