

**Product Name :**  
Universal Testing Machin Rotating Fatigue Machine

**Product Code :**  
NLAB-ENGINEERINGLB31003



**Description :**

Universal Testing Machin Rotating Fatigue Machine

**Technical Specification :**

Universal Testing Machin Rotating Fatigue Machine Lab Equipments Manufacturer, Suppliers & Exporters  
Loading accuracy as high + 1%. Speeds: Straining at variable speeds to suit wide range of materials. Facilities for tests: Motor-driven threaded columns for quick and convenient adjustment of lower cross head to facilitate rapid fixing of tear specimen. Autographic recorder: Simultaneous roll autographic recorder supplied as standard to enable study of the behaviors of materials. Ideal Dia.: High reading accuracy due to large size ideal design of dial. Large columns: Large effective clearance between columns enables testing of standard specimen as well as structures. Easy Changeability: Easy change from plain to threaded and screwed specimens. Simple and Safe: Simple to operate. Robust construction. Chrome plated metal components. Wide range of standard and special attachments / accessories available

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003