

**Product Name :**  
Universal Water Bath Serological

**Product Code :**  
NLAB-LBEMALAB230011



**Description :**

Universal Water Bath Serological

**Technical Specification :**

Matching the various IS standards are Shaker Water Bath made by a classic that reciprocates shaking motion application in the ever growing field related to agriculture, textile, medical, industry, research laboratories and hospitals.

Material: Stainless Steel

Shape: Rectangular

Ambient: + 5 degree C to 99 degree C

Motor: AC/DC Gear Motor

Accuracy: +/- 0.5 degree C

Revolutions: per minute 50 To 150

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003