

**Product Name :**  
Urinary Model Electronic

**Product Code :**  
NLAB-MEDICALMODELAB11007



**Description :**

Urinary Model Electronic

**Technical Specification :**

Leveraging upon our vast industry experience, we are offering a quality-assured range of Electronic Urinary Model.

It is manufactured in compliance with international standards using well-tested raw material.

This model is available in different sizes, diameter, length and thickness. We offer model in various technical specifications.

Clients can also avail our offered model in various customized options on their specific demands.

We offer this model in bulk quantities within promised time-frame to the clients.

Features:

Smooth finish

Washable

Perfect look

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003