

Product Name :
Urine Container

Product Code :
NLAB-EDUSCHOOLEQ133001



Description :

Urine Container

Technical Specification :

These containers can be used for collecting a wide range of samples of sputum.

Moulded in polypropylene these container have excellent contact clarity and the cap has built-in seal tdat makes it completely leakproof.

Urine Container Capacity 30ml.

Urine Container Capacity 60ml.

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com