

Product Name :
UV Chamber with UV Tube Light

Product Code :
NLABMEDICAL-M3190004

Description :

UV Chamber with UV Tube Light India, UV Chamber with UV Tube Light manufacturer India, suppliers, Medical UV Chamber with UV Tube Light India, Hospital UV Chamber with UV Tube Light India

Technical Specification :

The ultraviolet chamber we deal in is recommended for its various utility in all kinds of hospitals and laboratories.

These chambers are scientifically developed and designed by our highly engineers as per the international standard and norms.

They are accessible in high strength stainless steel body.

In addition to this our stainless steel UV chamber is backed with all modern and sophisticated features for maintaining the sterility of sterilized instruments in dental clinics.

We not only offer this chamber in customized order demanded by our clients but also at the most reasonable rates as against the general market.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003