

Product Name :
V Funnel Test Apparatus

Product Code :
NLAB -CONCRLAB120020



Description :

V Funnel Test Apparatus

Technical Specification :

The V Funnel apparatus is used to evaluate the flow time of freshly mixed selfcompacting concrete. The test is not suitable when the maximum size of the aggregate exceeds 22.4 mm.

The test set consists of a stainless steel funnel placed vertically on a supporting stand. The discharge orifice is equipped with a lid, which can be momentarily opened.

The I shape box apparatus is supplied complete with:

Filling Hopper

Base

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

