

Product Name :
V Funnel Test

Product Code :
NLAB-CONCTLAB240004



Description :

V Funnel Test

Technical Specification :

For determination of the relative falling time RM of mortar, completely made of Mild steel Powder Coated, with a mortar catch collar around the top.

Funnel w x d: 275 x 30 mm, height over all: 800 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003