

**Product Name :**  
Vaccine Carrier,Short Range

**Product Code :**  
NLABMEDICAL-M2320001



**Description :**

Vaccine Carrier,Short Range

**Technical Specification :**

Type: Vaccine Carrier- Short Range

Vaccine Storage Capacity As per Standard PQS load: 0.80 Litres

Gross Storage Capacity: 0.80 Litres

Weight-Fully Loaded: 2,18kg

Weight-Empty (with empty Ice Packs): 1,34kg

External Surface Material: HDPE (High Density Polyethylene)

Internal Lining Material: HIPS (High Impact Polyesterene)

Insulation Material: CFC Free Polyurethane

Insulation Thicknessv32-36mm

External Dimension (L X W XH): 25,00 x 18,00 x 21,00cm

Internal Dimension (L X W XH): 17,30 x 10,30 x 12,10cm

Vaccine Storage Dimensions: 17,30 x 10,30 x 04,50cm

Lid Type & Fixings: Removable

Number of Ice-packs Required: 2

Number of Ice-packs Supplied: 2

Ice-pack type: 0.4L compliant with PQS E0S/IP0I.1

Volume of Ice-pack: 0.4 Litres

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003