

**Product Name :**  
Vacuum Degassing Oven

**Product Code :**  
NLAB-ENGINEERINGLB31017



**Description :**

Vacuum Degassing Oven

**Technical Specification :**

The Vacuum Degassing Oven is used to precisely and accurately degas pressure-aged binder samples to meet standards.

The compact, table-top unit is constructed of stainless steel with a hinged lid to conserve space while allowing easy access to the stainless steel vacuum chamber. The oven holds up to 4 specimen containers and features a self-contained automatic vacuum system. The high pressure controller features a digital display indicating time, temperature and the current stage of each process, as well as illuminated pre-heat, degas and end/ stop buttons.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003