

Product Name :
Vacuum Pycnometer

Product Code :
NLAB-EDUCATIONTVET11016



Description :

Vacuum Pycnometer, Capacity 10 Liters

Technical Specification :

Dimensions Ø300 x h450 mm

Vibro-aerator

2-stage vacuum pump

Vacuum regulator

Water trap

Pipe

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

