

**Product Name :**  
Vacuum Tissue Processor

**Product Code :**  
NLAB-LABORATORYEM28035



**Description :**

Vacuum Tissue Processor

**Technical Specification :**

Designed for high capacity vacuum paraffin infiltration of tissue in routine clinical and research histopathology. A smart and user friendly software reduces user intervention that improves operating reliability. Superior process control and precise temperature regulation ensures optimum specimen quality. The smart start function allows one touch start of user-saved operations. It has a 200 Tissue cassettes capacity with a working cylinder temperature of 65°C and a volume of 5.3 litre. A large screen color LCD display, manual and automatic functioning mode, advanced safety features and a smart start function that allows one touch start of pre-saved programs. Three paraffin baths with a bath volume of 4.01 litre each, average paraffin melting time of 3 hrs. and operating temperature range of 15 to 35°C.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003