

Product Name :
Valve Seat Cutting Machine

Product Code :
NLAB-MACHINELABWK32006



Description :

Valve Seat Cutting Machine

Technical Specification :

Valve Seat Cutting Machine.

Formtools, small (20-30 mm dia) :

1. 15° x 45° x 75°-1 mm seal width
 2. 10° x 30° x 60°-1 mm seat width
 3. Counter boring tool
- 1 No

Formtools, -large (30-62 mm dia) :

1. 15° x 45° x 75°-1 mm seat width
2. 0° x 30° x 60°-1 mm seat width
3. Counter boring tool

STANDARD KIT :

The Standard Kit consists of the basic machine, power control unit, operating tools, carrying case, formtools and pilots as detailed.

WALL MOUNTING BRACKET :

Suitable for fixing all types of cylinderheads. The wall mounting bracket increases the efficiency of the seat cutting operation and reduces operator fatigue. This unit is essential for working on aluminium cylinder heads. (The bracket is supplied along With counter bore special mounting bolts.)

COUNTER BORING TOOL SETTING MICROMETER :

With the use of the optional Counter Boring Tool Setting Micrometer, the counter borecutting tools can be set to any specified diameter.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003