

Product Name :
Vamp Flex Tester

Product Code :
NLAB-TESTINGLABM94001



Description :

Vamp Flex Tester

Technical Specification :

Vamp Flex Tester. Use to find the resistance of a component of material to cut growth during repeated flexing, the offered flex tester is perfectly designed using

Features:

Easy installation

Rust resistance

Precisely designed

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003