

**Product Name :**  
Van Don Sampler

**Product Code :**  
NLAB-SCIENTIFICIND13009



**Description :**

Van Don Sampler

**Technical Specification :**

Water sampler ideal for water Column Sampling ;

Open water profiling ;

Sampling various stratification levels, the thermo cline or just above the bottom sediments;

Horizontal PVC includes carry case transparent PVC material 8.2L.

Can be used for shallow of deep applications;

Easy to clean and maintain ;

Can be used by one person in the field.

**Features & Benefits**

Can be used for shallow of deep applications

Easy to clean and maintain

Can be used by one person in the field

Sample can be seen by the user (acrylic versions only)

Applications

Water Column Sampling

Open water profiling

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003