


Product Name :
Variable Battery Internal Resistance Series

Product Code :
NLAB-SCHOOLABS13012



Description :

Variable Battery Internal Resistance Series

Technical Specification :

Internal resistance 5 – 30 Ω , Voltage 2 V. , Battery capacity 0.6 Ah , can prove Ohm's law

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003