

Product Name :
Variable Compression Equipment

Product Code :
NLAB-ENGINEERINGLB28002



Description :

Variable Compression Equipment

Technical Specification :

Description of Automotive Engine Test Bed Variable Compression Facility

The instrumentation provided is highly accurate, Includes an electric sensitive dynamometer, which allows To achieve fair and repeatable results. The engine can Be provided with a specific cylinder head fitted with Devices for the multi-channel indicator for Diagrams. Other cylinder heads for the gasoline engine Variable compression ratio are provided to allow To study the effects on the efficiency of the engine. Anything Instrumentation and controls are located on the Control panel made of melamine resin which comprises Multicolor synoptic table

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003