

Product Name :
Vee Bee Consistometer

Product Code :
NLAB -CONCRLAB120005



Description :

Vee Bee Consistometer

Technical Specification :

The instrument comprises of a slump cone with a hopper, specimen container, a transparent plastic plate attached to a graduated rod, mounted on a vibrating table, and a tamping rod 16mm dia. x 600mm long. The vibrating table has fixed amplitude and frequency of vibration, imparted by a motor drive mechanism

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003