

|   |   |
|---|---|
| <b>Product Name :</b><br>Velcro Cuff for Sphygmomanometer | <b>Product Code :</b><br>NLABMEDICAL-M3230005 |
|---|---|

|  |
|--|
|  |
|--|

**Description :**

Velcro Cuff for Sphygmomanometer India, Velcro Cuff for Sphygmomanometer manufacturer India, suppliers, Medical Velcro Cuff for Sphygmomanometer India, Hospital Velcro Cuff for Sphygmomanometer India

**Technical Specification :**

- Velcro cuff for sphygmomanometer, complete with rubber bladder with two tubes :For 1 - 3 years, size: 7 x 30.5cm..
- Velcro cuff for sphygmomanometer, complete with rubber bladder with two tubes :For new born, size: 5.7 x 30.5cm.
- Velcro cuff for sphygmomanometer, complete with rubber bladder with two tubes :For 3 - 7 years, size: 8.3 x 35cm.
- Velcro cuff for sphygmomanometer, complete with rubber bladder with two tubes :For 7 - 10 years, size: 9.5 x 38cm.
- Velcro cuff for sphygmomanometer, complete with rubber bladder with two tubes :For adult, sz: 14.5x 51.5cm.
- Velcro cuff for sphygmomanometer, complete with rubber bladder with two tubes :For obese, size: 17.5 x 63cm.
- Velcro cuff for sphygmomanometer, complete with rubber bladder with two tubes :For thigh, size: 21 x 83cm.

**Naugralabequipments**  
**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003