

Product Name : Velcro Cuff for Sphygmomanometer	Product Code : NLABMEDICAL-M3230005
Description : Velcro Cuff for Sphygmomanometer India, Velcro Cuff for Sphygmomanometer manufacturer India, suppliers, Medical Velcro Cuff for Sphygmomanometer India, Hospital Velcro Cuff for Sphygmomanometer India	
Technical Specification : Velcro cuff for sphygmomanometer, complete with rubber bladder with two tubes :For 1 - 3 years, size: 7 x 30.5cm.. Velcro cuff for sphygmomanometer, complete with rubber bladder with two tubes :For new born, size: 5.7 x 30.5cm. Velcro cuff for sphygmomanometer, complete with rubber bladder with two tubes :For 3 - 7 years, size: 8.3 x 35cm. Velcro cuff for sphygmomanometer, complete with rubber bladder with two tubes :For 7 - 10 years, size: 9.5 x 38cm. Velcro cuff for sphygmomanometer, complete with rubber bladder with two tubes :For adult, sz: 14.5x 51.5cm. Velcro cuff for sphygmomanometer, complete with rubber bladder with two tubes :For obese, size: 17.5 x 63cm. Velcro cuff for sphygmomanometer, complete with rubber bladder with two tubes :For thigh, size: 21 x 83cm.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003