

**Product Name :**  
Vernier Theodolite

**Product Code :**  
NLAB-ENGINEERINGLB85003



**Description :**

Vernier Theodolite

**Technical Specification :**

Vernier Theodolite least count 05/10/20 seconds as per customer requirements with four slow motion, erect image, 22x magnification, coated lens, along with powder coated/ anodized folding.

Aluminium tripod complete in plastic/ wooden carrying case.

Optical plummet arrangement is provided for centering.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003